

# QMRAspot: a tool for quantitative microbial risk assessment for drinking water

Manual QMRAspot version 2.0

RIVM report 2014-0020 J.F. Schijven et al.



# QMRAspot: a tool for quantitative microbial risk assessment for drinking water

Manual QMRAspot version 2.0

RIVM Report 2014-0020

# Colophon

#### © RIVM 2014

Parts of this publication may be reproduced, provided acknowledgement is given to: National Institute for Public Health and the Environment, along with the title and year of publication.

Jack Schijven, RIVM Saskia Rutjes, RIVM Patrick Smeets, KWR Watercycle Research Institute Peter Teunis, RIVM

Contact:
Saksia Rutjes
Z&O
Saskia.Rutjes@rivm.nl

This investigation has been performed by order and for the account of the Ministry of Infrastructure and the Environment, within the framework of the Knowledge Base Drinking Water (project number M/300005/01/CA)

This is a publication of:

National Institute for Public Health and the Environment P.O. Box 1 | 3720 BA Bilthoven The Netherlands www.rivm.nl/en

#### Abstract

# Manual for QMRAspot, a computational tool for the assessment of quantitative microbial risks of drinking water consumption.

The RIVM has developed a user-friendly computational tool (QMRAspot) to calculate the risk of becoming infected by pathogenic microorganisms in drinking water. This report is a manual in which it is explained how QMRAspot (Quantitative Microbial Risk Assessment from surface water to potable drinking water) can be used, and what the underlying models are. The tool is mainly aimed at pathogenic microorganisms in drinking water produced from surface water. The manual applies to the most recent tool version (2.0, September 2014).

The Dutch drinking water companies that produce drinking water from surface water and groundwater are obliged by law to demonstrate that less than one per then thousand persons per year acquire an infection by consumption of unboiled drinking water. QRMAspot was originally developed for the Netherlands, but can be applied worldwide by drinking water companies, researchers and policy makers.

This manual describes in detail how data for the risk assessments can be provided, how these data are analysed statistically, and how the risk assessments can be conducted in a consistent and transparent manner. In addition, two frequently occurring applications of QMRAspot are discussed. The first application demonstrates the effect on the infection risk of a high contamination at locations where the drinking water companies abstract surface water. The second application demonstrates to what extent one single sample, in which microorganisms were detected, contributes to the infection risk.

Keywords: Quantitative Microbial Risk Assessment; tool; drinking water; index pathogen

# Publiekssamenvatting

# Handleiding voor QMRAspot, een rekenprogramma voor het schatten van kwantitatieve microbiologische risico's van drinkwater.

Het RIVM heeft een gebruiksvriendelijk computerprogramma (QMRAspot) ontwikkeld dat de kans berekent op infecties door ziekteverwekkende microorganismen in drinkwater. Het onderhavige rapport is een handleiding, waarin wordt uitgelegd hoe QMRAspot (Quantitative Microbial Risk Assessment from surface water to potable drinking water) gebruikt kan worden en wat de onderliggende rekenmodellen zijn. Het model is voornamelijk bedoeld voor ziekteverwekkende micro-organismen in drinkwater dat uit oppervlaktewater wordt gewonnen. De handleiding is van toepassing op de meest recente modelversie (2.0, september 2014).

De Nederlandse drinkwaterbedrijven, die drinkwater produceren uit oppervlaktewater en grondwater, zijn wettelijk verplicht om aan te tonen dat minder dan één op tienduizend personen per jaar een infectie oploopt door de consumptie van ongekookt drinkwater. QRMAspot is oorspronkelijk ontwikkeld voor Nederland, maar kan wereldwijd worden toegepast door drinkwaterbedrijven, onderzoekers en beleidsmakers.

In deze handleiding wordt in detail beschreven hoe gegevens voor de risicoschattingen kunnen worden aangeleverd, hoe deze gegevens statistisch worden geanalyseerd en hoe de risicoschatting op consistente en transparante wijze kan worden uitgevoerd. Ook worden twee veel voorkomende toepassingen van QMRAspot besproken. De eerste toont het effect op het infectierisico van een hoge besmettingsgraad op locaties waar de drinkwaterbedrijven het oppervlaktewater onttrekken. De tweede toepassing demonstreert in welke mate één enkel monster waarin micro-organismen zijn aangetoond, bijdraagt aan het infectierisico.

Kernwoorden: Kwantitatieve microbiologische risicoschatting; tool; drinkwater; indexpathogeen

# Contents

1 1.1 1.2 1.3	$ \begin{array}{l} \textbf{Introduction} - \textbf{11} \\ \textbf{Drinking water legislation} - 11 \\ \textbf{General explanation of QMRA} - 11 \\ \textbf{Application of QMRA by drinking water companies} - 12 \\ \end{array} $
2 2.1 2.2 2.3 2.4 2.5	Excel spreadsheet with raw data $-$ 15 Raw data $-$ 15 Source water data $-$ 15 Recovery data $-$ 15 Treatment data $-$ 15 QMRAdata.xls $-$ 16
3.1 3.2 3.3 3.4 3.5 3.6	Run screen QMRAspot – 21 Run – 21 Select QMRA data spreadsheet – 22 Run QMRA – 23 Data analysis – 24 End of run – 25 Report – 26
4.1 4.2 4.3 4.4 4.5 4.6 4.7 4.8 4.9 4.10 4.11	Results: Tabs with index pathogen names $-$ 27 General $-$ 27 Data scheme $-$ 28 Source water $-$ 29 Source water with recovery correction $-$ 30 Treatment $z1 - 31$ Treatment $z2 - 32$ Treatment $z3 - 33$ Exposure $-$ 34 Infection risk $-$ 35 Summary table $-$ 36 Sensitivity analysis $-$ 40
5.1 5.2 5.3 5.4 5.5 5.6	Parameter settings – 41 General – 41 Source water – 42 Recovery – 43 Treatment – 44 Consumption – 45 Dose response – 46
6	Help — 47
<b>7</b> 7.1 7.2 7.3 7.4	Fitting distributions to the data $-$ 49 General $-$ 49 Source water concentration $-$ 49 Recovery, $R-50$ Treatment, $z-51$
8	Monte Carlo simulation — 53
9	Exposure and infection risk — 55

10	Dose response data – 57					
11	Example peak concentration in source water – 59					
12	Example treatment — 61					
12.1	General — 61					
12.2	Treatment z2, one positive sample in effluent — 61					
12.3	Treatment z2, only non-detects in effluent — 62					
12.4	Infection risk, one positive sample in effluent — 63					
12.5	Infection risk, only non-detects in effluent — 64					
12.6	Monte Carlo versus Point estimates — 65					
13	Contact – 67					
14	References — 69					

# Summary

QMRAspot is an interactive computational tool that has been designed for conducting Quantitative Microbial Risk Assessment (QMRA) for drinking water produced from surface water. This manual explains how to use QMRAspot as well as the underlying models in detail. QMRA results are explained on the basis of a reference data set. An example with the effect of a peak concentration in the source water and an example with treatment data consisting of only non-detects or only one positive sample in the effluent of the treatment step are shown. The manual applies to QMRAspot version 2.0 (1/9/2014).

In the Netherlands, drinking water companies are legislatively obligated to demonstrate compliance to not exceeding an infection risk of one per ten thousand persons per year from consumption of drinking water. This has to be demonstrated every four years for four index pathogens: Enterovirus, *Campylobacter, Cryptosporidium* and *Giardia*. QMRAspot has been designed to conduct QMRA for these index pathogens.

This report explains how QMRAspot assesses the infection risk and addresses possible applications. A comprehensive description of QMRAspot version 1.0 has been published by Schijven et al. (2011), which is fully cited in this manual, supplemented with explanation in more detail and examples, and updated to version 2.0 of QMRAspot. QMRAspot was originally developed in Mathematica 8.0.4 (Wolfram, Inc, Champaign IL, USA). The current QMRAspot version 2.0 has been updated to run with version 9.0.1 of Mathematica, Player Pro and CDF Player. QMRAspot version 2.0 has been extended to include a fifth optional pathogen, distribution parameter values for recovery can be set by the user, a dose response model can be selected, and dose response parameter values can be set by the user for the fifth pathogen. Raw data for OMRA to estimate the concentration of pathogens in the source water and to estimate removal of pathogens and/or indicator microorganisms by drinking water treatment must be stored in an Excel spreadsheet. These data can be read by QMRAspot with Mathematica or Player Pro. QMRAspot fits distributions to these data. In addition and/or instead, distribution parameter values can be set interactively in QMRAspot. A QMRA can thus be conducted without analyzing raw data as well. With the free CDF Player, data cannot be read from a spreadsheet, but QMRA can still be conducted by setting parameter values.

A QMRA report can be generated for each of the index pathogens, and this report can be saved as a Mathematica notebook and/or pdf file.

#### 1 Introduction

#### 1.1 Drinking water legislation

The World Health Organization (WHO) Guidelines for Drinking Water (WHO 2011) outline a preventive management framework for safe drinking water entailing health based targets, system assessment from source through treatment to the point of consumption, operational monitoring of the control measures in the drinking water production, management plans documenting the system assessment, and monitoring plans and a system of independent surveillance that verifies that the above are operating properly.

In line with the WHO Guidelines, the Dutch Drinking Water Act (2009) prescribes that tap water provided by the owner to consumers and other customers should not contain micro-organisms, parasites or substances to such numbers per volume or concentrations that these may comprise detrimental public health effects (Article 21). In the Dutch Drinking Water Act (2011), this demand is translated into the following quality requirements:

- 1) Absence of E. coli and enterococci in 100 ml of drinking water;
- 2) (Entero)viruses, *Cryptosporidium*, *Giardia* and *Campylobacter* should not exceed an infection risk of one infection per 10,000 individuals per year. In the Dutch Drinking Water Act (2011), no specific directions are given on how to perform this so-called Quantitative Microbial Risk Assessment (QMRA). To guide the owner of the provided tap water and the Inspectorate body on how to perform QMRA, the Inspectorate Guideline 5318 (Anonymous, 2005) was drafted in close consultation between the government (Environmental Inspectorate), the National Institute of Public Health and the Environment (RIVM), Bilthoven, the Netherlands, and the drinking water producers. The general principle of the Inspectorate guideline is to balance health protection and public funds.

#### 1.2 General explanation of QMRA

The required quantitative risk assessment is based on source water quality and the efficiency of the applied treatment. In addition, data are needed concerning tap water consumption and the dose-response relation of the specific pathogen and its host. Risk assessment for exposure to pathogenic microorganisms in drinking water was described by Teunis et al., 1997, Haas et al., 1999 and Haas and Eisenberg, 2001, the ILSI framework (Benford, 2001) and Medema et al. (2003). These publications and the development of risk based approaches by the WHO (2004) served as the basis for the specification of QMRA in Inspectorate Guideline 5318. It must be emphasized that risk assessment is an iterative process which is directed by practical and theoretical progress.

#### 1.3 Application of QMRA by drinking water companies

To demonstrate microbial safe drinking water, Dutch drinking water suppliers must conduct a QMRA at least every four years for the so-called index pathogens enterovirus, *Campylobacter*, *Cryptosporidium* and *Giardia*. Therefore, a three yearly monitoring program for index pathogens in the source water is prescribed in Inspectorate Guideline 5318, which may be condensed into one year to allow a higher monitoring frequency and hence provide more information about variability in pathogen concentrations. The Guideline also defines the monitoring frequency for source waters, which depends on the drinking water

production volume (Anonymous, 2005). In addition to regular monitoring, a number of incidental samples must be collected at moments when peak concentrations in pathogen counts are assumed to occur, for example due to heavy rainfall.

Because treatment efficiency is commonly highly location specific, treatment data should be collected at every production location. Any changes in the treatment process require a new estimation of the treatment efficiency, and thus new collection of data. Commonly, pathogen concentrations decrease below detection limits as a result of drinking water treatment. Therefore, indicator organisms are used to characterize treatment. They have similar properties as the index pathogens, and are assumed to be removed equally or less efficient by drinking water treatment. Appropriate indicator organisms occur in higher numbers and are easier to enumerate with higher recovery. Inspectorate Guideline 5318 (Anonymous, 2005) prescribes F-specific or somatic bacteriophages as the indicator organisms for determining removal efficiency by drinking water treatment of enterovirus. *Escherichia coli* is used as the indicator organism for *Campylobacter*, and spores of sulphite reducing clostridia (SSRC) are used for both *Cryptosporidium* and *Giardia*.

Because the efficiency of treatment varies in time, a sampling period should therefore be sufficiently extensive and frequent to be able to account for the most important sources of variation and should take changes in treatment into account. This also holds true for changes in the source water quality and changes in the scientific knowledge into the efficiency of treatment processes. Source water quality data need to be collected on a regular basis because of possible trends and year-to-year variations. Inspectorate Guideline 5318 (Anonymous, 2005) prescribes three year data collection for source water quality.

In order to automate the QMRA process, the interactive user-friendly computational tool, QMRAspot, was developed to conduct QMRA for drinking water produced from surface water. No extensive prior knowledge about QMRA modeling is required by the user, because QMRAspot provides the user with guidance on the quantity, type and format of raw data, and performs a statistical analysis of the raw data and then calculates a risk metric for drinking water consumption that can be compared with other production locations, a legislative standard, or an acceptable health based target. The uniform approach promotes proper collection and usage of raw data, warrants quality of the risk assessment, and improves efficiency, i.e., less time is required. QMRAspot facilitates QMRA for drinking water suppliers worldwide. The tool aids policymakers and other involved parties in formulating mitigation strategies, and prioritization and evaluation of effective preventive measures as integral part of water safety plans.

This report explains how QMRAspot assesses the infection risk and addresses typical applications. A comprehensive description of QMRAspot version 1.0 has been published by Schijven et al. (2011), which is fully cited in this manual, supplemented with explanation in more detail and examples, and updated to version 2.0 of QMRAspot. QMRAspot was originally developed in Mathematica 8.0.4 (Wolfram, Inc, Champaign IL, USA). The current QMRAspot version 2.0 has been updated to run with version 9.0.1 of Mathematica, Player Pro and CDF Player. QMRAspot version 2.0 has been extended to include a fifth optional pathogen, distribution parameter values for recovery can be set by the user, a dose response model can be selected, and dose response parameter values can be set by the user for the fifth pathogen. Raw data are stored in an Excel spreadsheet that can be read by QMRAspot with Mathematica or Player Pro. QMRAspot fits distributions to these data. In addition and/or instead, distribution parameter values can be set interactively in QMRAspot. A QMRA can thus be

conducted without analysing raw data as well. With the free CDF Player, data cannot be read from a spreadsheet, but QMRA can still be conducted by setting parameter values.

A QMRA report can be generated for each of the index pathogens and this report can be saved as a Mathematica notebook and/or pdf file.

# 2 Excel spreadsheet with raw data

#### 2.1 Raw data

For a QMRA, it is essential to collect quantitative microbial data as raw unprocessed data. Raw data on enumerated microorganisms in water are the counted numbers of the microorganisms as well as the corresponding investigated volume of the sample. Commonly, counts are numbers of plaqueforming units (pfu) for viruses and colony-forming units (cfu) for bacteria (Schets et al., 2008; Teunis et al., 2005a). Oocysts of *Cryptosporidium* and cysts of *Giardia* may be counted manually or automatically under a microscope using fluorescent dye (Schets et al., 2008). Raw presence/absence data are presence or absence of microorganisms in replicate dilutions of a water sample (De Roda Husman et al., 2009).

Obviously, the concentration of, for example, 1 pathogen particle in 1 millilitre of water is the same as a 100 particles in 100 millilitres of water, but 100 counted particles produce a more accurate concentration estimate than just one counted particle, hence it is essential to use counts and sample volumes as they were observed, and not concentrations (ratios of count/volume).

All raw data must include a sample date, which is needed to enable plotting time-series of the data. These plots may reveal variations, trends and extreme values, and may aid making selections of the data for QMRA in the spreadsheet.

# 2.2 Source water data

Source water data are raw data of index pathogens in the source water (e.g. Rutjes et al., 2009; Lodder et al., 2010). At many production locations, river water first passes a storage reservoir prior to further treatment. For enterovirus, river water may appropriately be designated as source water, because human contamination in a storage reservoir is not expected. For *Campylobacter*, *Cryptosporidium* and *Giardia*, the storage reservoir should be considered the starting point of the QMRA, because of contamination of the storage reservoir water with these pathogens from birds, wildlife, or runoff from agricultural land.

#### 2.3 Recovery data

In order to determine the recovery efficiencies of the detection method for index pathogens, ideally, each sample of source water, or a fraction used for analysis, is spiked with a sufficiently high number of, for example, a specific type of indicator organism. The spiked and recovered numbers can then be used to estimate the recovery efficiency. Therefore, raw recovery data consist of counts and samples sizes of the spiked and recovered microorganisms that are paired according to sampling date (Teunis et al. 1999; Rutjes et al., 2004; Schets et al., 2004). If the source water sample contains sufficiently high numbers of indicator organisms, like bacteriophages, to allow direct counting after plating, these indicator organisms can be used to estimate the recovery efficiency of a concentrating step needed for detection of the associated pathogen.

#### 2.4 Treatment data

Raw data of indicator organisms for a treatment process consist of counts and sample volumes of samples taken before and after treatment. If available, raw data of index pathogens can be used instead.

#### 2.5 QMRAdata.xls

The tool reads the raw data from a standard Microsoft Excel spreadsheet file, here, for convenience, called QMRAdata.xls. It contains three sheets: SCHEME, RAW DATA and HELP. The SCHEME sheet provides a description of the drinking water production location and defines a table with column headers that is used by the tool to make the appropriate data selections. Through the SCHEME sheet, the user has control over what data should be used for QMRA. Obviously, the RAW DATA sheet contains all raw data and the HELP sheet provides background information on how to fill the RAW DATA sheet with raw data in the required format.



Figure 2.1 SCHEME sheet of QMRAdata.xls

In the SCHEME sheet, the effluent data of a treatment step can be the influent data for the next treatment step. The data for the next treatment step may also be from other types of microorganisms, or from a different location in the drinking water utility, or from pilot plant experiments.

It is easy to modify the SCHEME sheet. For example, two treatment steps (designated z1, z2) may be combined into a single one by using the influent data of the first of the two and the effluent data of the second of the two treatment steps; in this example, z1+z2 with the influent of z1 and the effluent of z2. This can be done for each index pathogen independently.

The SCHEME sheet is a form in which the names of the drinking water company, the production location, and the names of the treatment steps can be given. Also, the names or labels for the influent and effluent of each treatment step need to be given.

The names of the source water and of the influents and effluents need to be exactly the same names as in column A Item in RAW DATA. QMRAspot uses these names for selecting the associated data.

The source of the data for the treatment steps can be selected from the following list: Plant scale, pilot plant scale, laboratory scale, and literature. Location-specific plant scale data are generally preferred, followed by pilot plant scale data, and if these are not available, data from laboratory experiments. In other words, location specific data are recommended. If the use of data from other locations is desired, applicability should be verified by comparison of treatment conditions. References to data from literature should be listed in the accompanying QMRA reports. How data from literature can be entered is described in chapter 5.

The SCHEME includes the four index pathogens Enterovirus, *Campylobacter*, *Cryptosporidium* and *Giardia*. By default, the corresponding indicator organisms for recovery are F-specific RNA bacteriophages or somatic coliphages for Enterovirus, *E. coli* for *Campylobacter* and SSRC (Spores of Sulphite Reducing Clostridia) for *Cryptosporidium* and *Giardia*, but other (arbitrary names of) microorganisms are also possible.

The names of the index pathogens and indicator microorganisms for recovery and treatment in RAW DATA column C Microorganism need to be exactly the same as given in SCHEME, so that QMRAspot can select the associated data.

In version 2.0 of QMRAspot, a fifth optional pathogen can be added, for which dose-response data can be entered interactively (see section 5.6). This can be any pathogen. QMRAspot version 2.0 retains the ability to read older spreadsheets that contain data for the four default index pathogens only.

	J 47 - CH -		A STATE OF	all and a second	1236									
- 1	File Home		e Layout Formu	las Data Review	v View F	Maces	Acrob	at						
	X10	<b>→</b> (*)	fx											
d		A	В	С	D	E	F	G	Н	1	J	K	L	N
1	Label	~	Sampling cod *	Microorganism	▼ Date ▼	Cc ₹			R1 *	MPN1 -	V2(liter) ▼	R2 *	MPN2 *	V3(lite
2	Source water			Enterovirus	1/7/2008	1	546							
3	Source water			Enterovirus	1/14/2008	0	508.3							
4	Source water			Enterovirus	1/21/2008	0	556							
5	Source water			Enterovirus	1/28/2008	1	520.2							
6	Source water			Enterovirus	2/6/2008	0	496.6							
7	Source water			Enterovirus	2/11/2008	0	504.6							
8	Source water			Enterovirus	2/18/2008	0	556.8							
9	Source water			Enterovirus	2/25/2008	0	495.1							
10	Source water			Enterovirus	3/3/2008	0	521.6							
11	Source water			Enterovirus	3/10/2008	1	530.5							
12	Source water			Enterovirus	3/17/2008	0	522.8							
13	Source water			Enterovirus	3/26/2008	0	502.3							
14	Source water			Enterovirus	3/31/2008	0	525.2							
15	Source water			Enterovirus	4/7/2008	1	512.7							
16	Source water			Enterovirus	4/14/2008	0	532.9							
17	Source water			Enterovirus	4/23/2008	0	539.2							
18	Source water			Enterovirus	5/14/2008	0	520.7							
19	Source water			Enterovirus	5/28/2008	0	511.5							
20	Source water			Enterovirus	6/9/2008	0	499.5							
21	Source water			Enterovirus	6/23/2008	0	504.2							
22	Source water			Enterovirus	7/7/2008	0	500.4							
23	Source water			Enterovirus	7/21/2008	0	503.7							
24	Source water			Enterovirus	8/4/2008	0	441.8							
25	Source water			Enterovirus	8/8/2008	0	552							
26	Source water			Enterovirus	9/1/2008	0	538.7							
27	Source water			Enterovirus	9/15/2008	0	544.5							
28	Source water			Enterovirus	9/29/2008	0	515.5							
29	Source water			Enterovirus	10/13/2008	0	515							
30	Source water			Enterovirus	10/27/2008	0	530.9							
31	Source water			Enterovirus	11/10/2008	0	560.6							
32	Source water			Enterovirus	11/17/2008	0	481.4							
33	Source water			Enterovirus	11/24/2008	0	526.6							
34	Source water			Enterovirus	12/1/2008	0	525.6							
35	Source water			Enterovirus	12/8/2008	0	528.4							
36	Source water			Enterovirus	12/15/2008	1	473							
37	Source water			Campylobacter	1/8/2007	40	1							
38	Source water			Campylobacter	1/15/2007	110	1							
39	Source water			Campylobacter	1/22/2007	9	1							
40	Source water			Campylobacter	1/29/2007	240	1							
41	Source water			Campylobacter	2/5/2007	240	1							
42	Source water			Campylobacter	2/12/2007	40	1							
43				Campylobacter	2/19/2007	50	1							
44				Campylobacter	2/26/2007	23	1							

Figure 2.2 RAW DATA sheet of QMRAdata.xls

In the RAW DATA sheet, every row is a full record of raw data. This simple design allows for automated filling from a Laboratory Information Management System (LIMS), for example, as records are stored in the form of Comma Separated Values (CSV).

- Column A Label: The name of source water, "Spike", "Recovered", treatment influent, or treatment effluent.
- Column B Sampling code: Specific code. Not used in QMRA calculation, but included in RAW DATA to allow for more detailed reference. Different codes may be included for different sampling points of the same effluent. In the OMRA these are combined.
- Column C Microorganism: Index pathogens and indicator microorganisms. Exactly the same names as given in SCHEME.
- Column D Date: Date format: DD-MM-YYYY.
- Column E Count: Only positive integers are allowed. Max counts per plate
  according to standard method. Viral counts are counts of enterovirus or
  bacteriophage plaques. Bacterial counts are colony counts in detection
  methods for bacteria using membrane filtration. Protozoal counts are
  microscopic (automated) counts of fluorescent labelled (oo)cysts.
- Column F Sample size: Dimension liter. The original equivalent volume of the sample that the detected microorganisms were counted in.
   Example 1: a 10-litre sample was collected, concentrated to 100 ml, 5 ml was plated for counting, then the sample size is 10/20=0.5 litre.
   Example 2: 98 colonies were counted on a plate with 1 ml sample and 11 colonies were counted in the ten-fold dilution. Count is 109 and sample size is 0.0011 litres.

• Columns G-O: MPN data, Most Probable Number data. May be available instead of counts and sample volumes.

V1, V2, V3, V4, V5 are the sample volumes of five dilution steps. Dimension is litre.

R1, R2, R3, R4, R5 are the numbers of replicates for each dilution step. Only whole numbers allowed.

MPN1, MPN2, MPN3, MPN4, MPN5 are the numbers of positive samples of a particular analysed volume. Only whole numbers from 0 - the number of replicates are allowed.

#### Example

V1(liter)	R1	MPN1	V2(liter	R2	MPN2	V3(liter)	R3	MPN3
0.001	3	3	0.0001	3	2	0.00001	6	4

In this example, the volume of dilution step 1 is 0.001 litres and it is replicated three times. Of those replicates, all three showed positive detection.

Recovery: RAW DATA of "Spiked" and "Recovered" microorganisms. The label is "Spike" or "Recovered". The data are counts and sample sizes.

Treatment: In SCHEME, the name of the treatment step can be given. This may also be a combination of treatment steps.

# 3 Run screen QMRAspot

#### 3.1 Run




Figure 3.1 Run screen of QMRAspot at start up

At start up, the Run screen of QMRAspot has no file selected; the index pathogens Enterovirus, *Campylobacter*, *Cryptosporidium* and *Giardia* are all selected. Obviously, QMRAspot only conducts QMRA for the index pathogens that are selected.

Pressing the "Reset" button resets QMRAspot to this state. Resetting is disabled while QMRAspot runs.

Should a user want to stop a run of QMRAspot by aborting an evaluation, it is recommended to restart with a new kernel or reloading of QMRAspot.

# 3.2 Select QMRA data spreadsheet




Figure 3.2 Run screen of QMRAspot after selecting spreadsheet with QMRA data

The "Select file" button allows the selection of a QMRAdata spreadsheet. In this case QMRAdataWRES.xls has been selected. This spreadsheet contains a combination of data from several drinking water companies (Schijven et al., 2011). This particular data set can be used as the reference data set. It serves as a template for other QMRA data sets and can be used to compare QMRA outcomes calculated by different versions of QMRAspot to identify possible differences. Together with QMRAspot, each new user should have a copy of the reference data set.

# 3.3 Run QMRA




Figure 3.3 Run screen of QMRAspot at fitting distributions

The QMRAspot logo is a button. Pressing it starts the QMRA run. On the start-up screen, a progress bar allows the user to monitor that after a few seconds, QMRAspot has read all the data on source water concentration, recovery and treatment from the spreadsheet and is fitting the appropriate distribution to the data of each QMRA step.

# 3.4 Data analysis




Figure 3.4 Run screen of QMRAspot almost finished

After a short while (in this case 1 minute and 6 seconds), QMRAspot has finished fitting distributions, plotting time series of the data read from the spreadsheet, drawing samples from the fitted distributions (for Monte Carlo calculations), and overall treatment, exposure and infection risks have been calculated.

#### 3.5 End of run




Figure 3.5 Run screen of QMRAspot at end of run

After finishing all the calculations, the main QMRAspot screen shows boxwhisker plots of estimated infection risks per person per year for the selected index pathogens. The y-axis shows the infection risk per person per year on a  $\log_{10}$  scale. The whiskers of the box-whisker plots are the 5- and 95-percentile values of the infection risk, and the box entails the 25- and 75- percentiles. The (arithmetic) mean infection risk is indicated with a blue line. The calculated risk is presented in relation to the Dutch health target of  $10^{-4}$  per person per year. If the mean AND the 95-percentile lie below  $10^{-4}$  per person per year (indicated as a black horizontal line) then the box is coloured green to denote compliance with the health based target. The box is orange if the mean OR the 95-percentile exceeds the target. If the box is red, then the mean infection risk AND the 95-percentile exceed the target. If the 95-percentile of the infection risk is lower than  $10^{-9}$  or higher than 0.9, a message appears with this information instead of a box-whisker plot.

Note that a proper risk assessment outcome needs to include a risk value and its probability. In this QMRA, the mean infection risk does not fulfil that requirement. The probability corresponding to any infection risk depends on the shape of its distribution. It is most straightforward to use the 95-percentile value as the target. This choice implies that per person per year, a predicted infection risk of  $10^{-4}$  is not exceeded with 95% confidence.

# 3.6 Report




Figure 3.6 Run screen of QMRAspot at end of run with QMRA report for enterovirus in a new window.

The box-whisker plots on the run screen that appear at the end of the run are buttons that when pressed open a window with a QMRA report of the selected pathogen. The report, or selections in the report can be printed and saved as pdf.

Each report contains a general table with reference to the tool version, run date, QMRA data Excel file and which index pathogen. A drinking water utility table gives a summary description of the drinking water production location. The summary table and all histograms are included.

# 4 Results: Tabs with index pathogen names

#### 4.1 General

Under the tabs with the index pathogen names the result of all the steps in the QMRA show up in detail as soon as they are available.

In the different stages of the QMRA, histograms of ten thousand Monte Carlo simulated data are presented; see section 8 for a brief explanation of Monte Carlo simulation. Ten thousand Monte Carlo samples have empirically been found to be a sufficient number of samples for a stable outcome (Teunis and Havelaar, 1999).

Figure 4.1 shows a histogram of Monte Carlo data as an example. These data are Monte Carlo samples that were drawn from a Beta distribution with parameters  $\alpha$ =0.24 and  $\beta$ =970, and represent ten thousand fractions of microorganisms passing a treatment step. The left histogram shows the probability distribution of the fractions. As can be seen, most of the fractions are very small values; therefore, a  $\log_{10}$  transformation of the fractions is appropriate. The left histogram shows the probability distribution of the  $\log_{10}$ -transformed fractions from the 1- percentile to the 99-percentile. To the right of this histogram is a table presenting information on the type of distribution, distribution parameter values and some obvious statistics. Note that this distribution spans a large range. Values of -8  $\log_{10}$  or lower are in fact of numbers near zero, all with a small probability.




Figure 4.1 Histogram with Monte Carlo data

# 4.2 Data scheme




Figure 4.2 Enterovirus-Data scheme

Under Enterovirus-Data scheme, a summary table is shown of all data that were read from the spreadsheet. It shows the names of the source water and treatment steps, and from what indicator organism treatment data were available. The rightmost two columns show the numbers of samples available for each QMRA step.

#### 4.3 Source water




Figure 4.3 Enterovirus-Source water, no recovery data

Under Enterovirus-Source water, the source water data are presented. The top left plot shows all individual pathogen concentrations in the source water as a time series. In this particular example, in five samples (#pos) of 35 samples (#tot), enterovirus was detected. The blue line is the mean concentration of the fitted Gamma distribution for the pathogen concentration. The area shaded in light blue encompasses the area below the 95-percentile of the fitted Gamma distribution data. In five cases, corresponding to the positive samples, the estimated concentrations appear to exceed the 95-percentile value. By definition, these values may be considered to be peak concentrations for which it is recommended investigating what environmental conditions may explain the occurrence of these high values, if not already known.  $\Sigma N$  and  $\Sigma V$  represent the total count and sample volume.

The bottom left panel shows a histogram of ten thousand Monte Carlo values (see also section 8) that were drawn from the fitted Gamma distribution that describes the variability of the pathogen concentrations in the source water, and a table with parameters r and  $\lambda$ , and some statistics of that distribution. The bottom right panel shows the distribution of the pathogen concentration after correction for recovery. In this example, it is the same histogram as in the left hand panel. In this example, no recovery data were available; hence, it was assumed that recovery was equal to one. So far, lack of specific recovery data is common, unfortunately, and it should be realized as well as reported that in this case the infection risk may be underestimated.

#### Run Parameter settings Enterovirus Campylobacter Cryptosporidium Giardia Help Data scheme | Source water | z1 | z2 | z3 | z4 | z5 | z6 | Exposure | Infection risk | Summary table | Sensitivity analysis Source concentration 0.06 0.05 0.0023 0.06 0.04 95% 0.0085 N/L 0.03 P 0.04 Mean 5% 0.033 0.02 0.02 0.01 0.00 Log<sub>10</sub>N/L Source concentration 0.12 ≤10<sup>-9</sup> ≤10<sup>-9</sup> 0.0023 5% ≤10<sup>-9</sup> ≤10<sup>-9</sup> -9.5 5% -10. ----III Log<sub>10</sub>N/L Log<sub>10</sub>N/L

# 4.4 Source water with recovery correction

Figure 4.4 Cryptosporidium-Source water with recovery data

This screen shows the source water data for *Cryptosporidium* including recovery data. The recovery data encompass paired data of spiked and recovered numbers of microorganisms to which a Beta distribution with parameters  $\alpha$  and  $\beta$  is fitted. A Beta distribution describes a variability of a fraction, *i.e.* a number between 0 and 1. The top right panel shows the fitted distribution as a histogram of ten thousand Monte Carlo samples. The ten thousand random samples from the Gamma distributed pathogen concentration in the source water (bottom left panel) is divided by the ten thousand random Beta distributed recovery fractions, to produce ten thousand corrected source water concentrations, shown in the bottom right panel, as a histogram and a table with a few statistics; see section 8 for a brief explanation of Monte Carlo simulation. The distribution of the corrected concentration has shifted approximately one log<sub>10</sub> to the right and is a bit wider than the distribution of the source concentration.

#### 4.5 Treatment z1




Figure 4.5 Enterovirus-z1, the first treatment step "Step 1"

Enterovirus-z1 shows influent and effluent somatic coliphage concentrations of the first drinking water treatment step z1 as a time series, with the 95% region shaded in light blue and the average concentration as a blue line, similar to the graphs of the pathogen concentrations in the source water. The time series plot on the bottom left shows both influent and effluent concentrations on a logarithmic scale. In that combined plot, non-detects cannot be included. The purpose of this plot is to allow the user to scrutinize possible variation of the log removal of microorganisms by treatment. This example shows that both influent and effluent concentrations vary by season, and that removal efficiency varies as well.

Unlike recovery data, treatment data are not treated as paired, even when influent and effluent were sampled on the same day, *i.e.* the microorganisms that are counted in the an effluent sample may not have originated from the same influent as the microorganisms that were counted in the influent sample of the same date. Also, often the numbers of samples taken from influent and effluent are not even the same. Obviously, influent and effluent data must be of the same time period. The bottom right panel shows the variation in removal, as a histogram of ten thousand beta distributed fractions of microorganisms that passed this treatment stage, and some statistics (see section 7.4). Note that reduction is expressed on a  $\log_{10}$  scale: in this example, on a  $\log_{10}$  scale, the distribution is slightly left skewed.

Note that in the tab "z1" is designated as a general reference for the first treatment step, and that "Step 1" is the name of treatment step z1 as read from the spreadsheet.

#### 4.6 Treatment z2




Figure 4.6 Enterovirus-z2, the second treatment step "Step 2"

Enterovirus-z2 shows influent and effluent somatic coliphage concentrations of the second drinking water treatment step z2 as a time series (top panels and bottom left). In this case, half of the effluent data are non-detects. The beta distributed fractions (bottom right panel) are left skewed.

# 4.7 Treatment z3




Figure 4.7 Cryptosporidium-z3, the third treatment step "Step 3"

*Cryptosporidium*-z3 shows influent and effluent SSRC concentrations of the third drinking water treatment step z3 as a time series. In this case, 538 of 582 effluent data are non-detects. The beta distributed fractions are strongly left skewed.

#### 4.8 Exposure




Figure 4.8 Enterovirus-Exposure

The top left panel of the exposure screen shows the combined treatment effect as a histogram of log reductions, and a table with statistics. For each of the treatment steps there are ten thousand randomly sampled fractions. The first fractions of each of the treatment steps are multiplied with each other, and so are the second, third, etcetera. This produces the total treatment as shown in the top left histogram. By multiplying the corrected pathogen concentration data with the total treatment data, the drinking water concentration data are obtained, shown in the bottom left panel.

The top right panel shows tap water intake as the (lognormal) distribution of consumption of un-boiled drinking water per person per day in the exposed population. Multiplication of drinking water concentration data with these consumed volumes produces the exposure, or the dose, as the numbers of ingested pathogens per person per day as a distribution in the bottom right panel.

## 4.9 Infection risk




Figure 4.9 Enterovirus-Infection risk

The infection risk screen shows the infection risk per person per day (top left panel) calculated by applying the dose response model to the exposure estimates. The yearly infection risk (risk of infection per person per year) is calculated by repeated sampling from the daily infection risks as explained in chapter 9 (top right panel). The infection risk per person per year is also shown in a box-whisker plot (bottom centre).

#### 4.10 Summary table

The Summary Table of QMRA results is given for all four index pathogens using the reference data set. These tables show mean, 5-percentile, median and 95-percentile values for each QMRAstep including their  $\log_{10}$  values. For comparison, point estimates are included in the rightmost column. The point estimate for the pathogen concentration in the source water is calculated from the total counted number and total volume of all samples. The point estimate is the maximum likelihood value of Poisson-distributed counts, assuming a fixed (constant) concentration. Such weighted average concentrations are also calculated for the indicator organisms in the influent and effluent of each treatment step. Their ratios represent a point estimate for the fraction of microorganisms that pass the corresponding treatment stage(s). Comparison of the results from the Monte Carlo simulation and point estimates demonstrates that the point estimate of exposure may be higher (0.5  $\log_{10}$ ) or lower (1  $\log_{10}$ ) than the mean exposure calculated from the Monte Carlo simulation.

Note that the use of a point estimate ignores information about its probability. In Monte Carlo simulation, any sampled number occurs with its own probability as demonstrated by the quantiles below.




Figure 4.10 Enterovirus-Summary table

The Enterovirus summary table shows that the point estimate of the exposure is  $0.3 \log_{10}$ , i.e. a factor of two lower than the mean exposure. This difference mainly resides from treatment step 2 (section 4.6).




Figure 4.11 Campylobacter-Summary table

In the *Campylobacter* summary table, the point estimate of the pathogen concentration in the source water is  $0.5 \log_{10}$  lower than the mean of the distribution. Point estimates of all treatments steps suggest more efficient removal than the Monte Carlo simulation, leading to almost  $1 \log_{10}$  lower exposure according to the point estimates than with the Monte Carlo sampling approach. Note that the mean exposure as calculated by Monte Carlo simulation exceeds the 95-percentile because of the skewness of the distribution (Figure 4.11).




Figure 4.12 Campylobacter-Exposure histogram with mean value larger than the 95-percentile

Data scheme   Source water   z1			z2	z3	z3 z4 z5 z		6 Exposu	ıre Infection risk		Sur	nmary table	Sensitivity analysi			
1									aidiaa DIV	M D:IAL		- 150530	11.4		
	Cryptosporidium RIVM Bilthoven 150 OMRA-step Dimension Mean 5% M									Media		95%	Point Estimate		
					+		N/L	-	0.0022		:10 <sup>-9</sup> ≤10 <sup>-9</sup>			0.01	0.0023
	Source concentration Recovery					N/∟ ction		0.0022	0.04		0.1		0.01	0.0023	
								<10 <sup>-9</sup>	<10 <sup>-9</sup>						
	Corrected source concentration			n			0.018	≤10				0.07	0.014		
	Step 1					fraction		0.35					0.69	0.38	
	Step 2					fraction		0.0042					0.012	0.0045	
	Step 3				fraction		0.031	0.000				0.21	0.055		
	Total treatment					fraction			0.000046	≤10 <sup>-9</sup>		$2.5 \times 1$		0.00022	0.000096
	Drinking water concentration				N/L		$8.1 \times 10^{-7}$	≤10	≤10 <sup>-9</sup> ≤10 <sup>-9</sup>		9	6.7×10 <sup>-8</sup>	1.3×10 <sup>-6</sup>		
	Consumption			L	L/person/day		0.28	0.02	26	0.16		0.9	0.28		
	Exposure			N	N/person/day		$1.7 \times 10^{-7}$	≤10	9	≤10	9	1.×10 <sup>-8</sup>	$3.7 \times 10^{-7}$		
	Infection risk				/pers	on/da	av I	4.2×10 <sup>-8</sup>	0		0.		2.2×10 <sup>-9</sup>		
	Infection risk			- 1 '		on/ye		0.000015	1.×1	0-6	7.1×1	0-6	0.000072		
					$\equiv$	Dimension									
	QMRA-step						n							Log <sub>10</sub> Point Estimate	
	Source concentration			N/L			-2.7	-42		-10		-2.	-2.6		
	Recovery					ction	۱ ا	-0.77	-1.3		-0.8		-0.48	-0.78	
	Corrected source concentration			n		N/L		-1.7	-41		-9.6		-1.2	-1.9	
	Step 1					ction	- 1	-0.46	-1.		-0.4		-0.16	-0.42	
	Step 2					ction	- 1	-2.4	-3.		-2.5		-1.9	-2.3	
	Step 3 Total treatment				fraction			-1.5	-17		-4.4		-0.69	-1.3	
								۱ ا	-4.3	-20		-7.6		-3.7	-4.
	Drinking water concentration				N/L			-6.1	-52		-19		-7.2	-5.9	
	Consumption				L/person/day N/person/day		-0.55			-0.8		-0.046	-0.55		
	Exposure							-6.8	-52		-20	•	-8.	-6.4	
	Infection risk					on/da		-7.4 -4.8	-6		-5.2		-8.7 -4.1		
	Infection risk					/person/year							-4.1		

Figure 4.13 Cryptosporidium-Summary table

In the *Cryptosporidium* summary table, the point estimate of the exposure is  $0.4 \log_{10}$  higher than the mean from the Monte Carlo simulation, mainly because of treatment step z3. Note that, here again, the mean exposure as calculated by Monte Carlo simulation is higher than the 95-percentile because of the skewness of the exposure distribution (Figure 4.13).




Figure 4.14 Cryptosporidium-Exposure histogram with mean value larger than the 95-percentile




Figure 4.15 Giardia-Summary table

In the *Giardia* summary table, the same indicator organism is applied for treatment efficiency as for *Cryptosporidium* and, therefore, overestimation by the point estimates is also similar to that of *Cryptosporidium*.

## 4.11 Sensitivity analysis




Figure 4.16 Enterovirus-Sensitivity analysis

To illustrate the importance of variable factors in QMRA, a simple sensitivity analysis is conducted in which the contribution of each step to the variance of the exposure is calculated. To that end, all Monte Carlo estimates are log transformed, and for each step in QMRA, the variance is calculated and divided by the variance of the exposure or dose. The sensitivity graph in the above example shows these contributions, sorted by their magnitude. It appears that treatment Step 2 accounts for the main contribution to the variance in exposure. On the basis of this sensitivity analysis, the drinking water company may decide to reduce variation in operational conditions of treatment Step 2, as part of their water safety plans, thereby increasing reliability and enabling better risk prediction.

## 5 Parameter settings

#### 5.1 General

Instead of and/or in addition to reading raw data from a spreadsheet, it is also possible to simulate pathogen concentrations in source water as well as removal of microorganism by treatment by entering characteristic distribution parameter values and then conducting a QMRA. If QMRAspot runs with the CDF Player, entering characteristic parameter values is the only way of conducting the QMRA, because in that case spreadsheet data cannot be imported. Although raw data are the preferred basis for QMRA, the option of QMRAspot for entering distribution parameters directly enhances its versatility. QMRA can thus also be conducted, if raw data are lacking, or it can be used to answer a variety of whatif questions, or do scenario studies.

While a QMRA is running, parameter settings cannot be changed. Once a spreadsheet with QMRA data is selected, names of source water, microorganisms and treatment steps are taken from the spreadsheet and cannot be set or altered under parameter settings.

#### 5.2 Source water




Figure 5.1 Parameter settings-Source water

Under Parameter settings-Source water, for each selected index pathogen, it is possible to have the tool read raw data, or to enter a concentration mean  $\mu$  and 95-percentile (can be approximated with a maximum value), or to use parameter values for r (shape parameter) and  $\lambda$  (scale parameter) of a Gamma distribution. Commonly, one may have at least some notion of a mean and maximum source water concentration. When these are given, the tool will directly estimate the corresponding Gamma distribution parameters. Note that it is a property of the Gamma distribution that its 95-percentile value cannot exceed the mean by more than 5.8 times, which nevertheless represents a very wide distribution. When the Gamma distribution parameters are given, QMRAspot calculates the corresponding mean and 95-percentile. Suppose a few samples have been analysed for one of the index pathogens. A proper estimate of the mean concentration is the weighted average concentration as calculated from the total counts and sample size of all samples. The 95-percentile is usually the highest concentration.

Estimates for r and  $\lambda$ , or mean and 95-percentile may also be taken from literature, or even be fictive values used for scenario calculations. If no spreadsheet with QMRA-data has been selected, then names for the source waters can be entered by the user.

The raw data option is not available when QMRAspot runs with the CDF Player.

## 5.3 Recovery




Figure 5.2 Parameter settings-Recovery

Under Parameter settings-Recovery, for each index pathogen, it is possible to set the tool to reading raw data, or using mean recovery R and parameter a of the Beta distribution, or using parameter values for a and  $\beta$  of the Beta distribution. R, a, and  $\beta$  need to be larger than zero. If R=1, or input of R, or a, or  $\beta$  is not a positive number, then R will be set to one. If values for a and  $\beta$  are available, for example from literature, their values can be entered directly. If one has some notion of mean and 95-percentile, then first the mean value needs to be entered directly, and then a value for parameter a needs to be entered such that the desired 95-percentile value is obtained.

If no spreadsheet with QMRA-data has been selected, then names for the indicator organisms can be entered by the user.

The raw data option is not available when QMRAspot runs with the CDF Player.

#### 5.4 Treatment




Figure 5.3 Parameter settings-Treatment-Enterovirus

Although QMRA based on location-specific raw data for treatment at full scale is strongly preferred, the tool provides the option for including distribution parameters values of fraction z of the microorganisms that were able to pass treatment instead of raw data. This option is not included to move away from collecting raw data, but often location-specific data at plant scale are not available because indicator organism levels were (expected to be) below detection limits. This often occurs with very efficient treatment steps and/or at the end of the production chain. In those cases, one has to rely on (literature) data from pilot plant experiments that mimic full scale conditions, or on data from laboratory scale experiments, or use treatment data that were collected at other plants operating under similar conditions. In all these cases, the applicability of the data to the location specific conditions needs to be verified. It is also possible to use models for treatment processes to predict removal values or to provide Beta distribution parameter values.

The option of including a treatment step by means of its distribution parameters can also be used to determine the required additional treatment if a drinking water location exceeds the health based targets. This option allows for scenario studies, and therefore, greatly increases the versatility of QMRAspot. Under Parameter settings - treatment, for each index pathogen, a mean log<sub>10</sub> value for log removal by treatment can be set instead of using raw data. In CDF Player, this is the only option. By setting the value of Beta distribution parameter a, the shape of the Beta distribution can be set. The associated values of the 95-percentile and Beta distribution parameter  $\beta$  are calculated and given. The user can experiment with these settings to achieve their desired target values. If input of  $\log_{10}z$ , or a, or  $\beta$  is not a positive number, then  $\log_{10}z$ will be set to zero, implying no treatment; see also the previous section on how to enter the parameter values. It is also possible to switch off a treatment step. If no spreadsheet with QMRA-data has been selected, then names for the treatment steps and indicator organisms can be entered by the user. The raw data option is not available when QMRAspot runs with the CDF Player.

## 5.5 Consumption




Figure 5.4 Parameter settings-Consumption of unboiled drinking water

QMRAspot offers four alternatives for consumption of unboiled drinking water per person per day. Except for the third option, these are all lognormal distributions, with parameters defined by various studies. Parameters  $\mu\text{=-}1.85779$  and  $\sigma\text{=}1.07487$  are for the Netherlands, corresponding to a mean of 0.27 litres per person per day (Teunis and Havelaar, 1999), a lognormal distribution with parameters  $\mu\text{=-}0.03598$  and  $\sigma\text{=}0.77218$  for the USA, corresponding to a mean of 1.3 litres per person per day (USEPA, 2006), a fixed volume of 2 litres per person per day (WHO, 2011) and, finally, the possibility of putting in the parameter values for  $\mu$  and  $\sigma$  for any other lognormal distribution of drinking water consumption, if available for another country or for a specific subpopulation. Consumption data may differ between countries and also between subpopulations; climate may also play a role. More data and a discussion about such variability, is given in USEPA (2006), by Westrell et al. (2006a) and WHO (2008).

## 5.6 Dose response




Figure 5.5 Parameter settings-Dose response

Four dose response models are available and shown in the table.

- 1. Exposure probability in the case no dose response data are available.
- 2. The exponential model, with parameter r. If r=1 then it is assumed that probability of infection equals exposure probability (option 1.).
- 3. The exact Beta Poisson model using the hypergeometric function  ${}_1F_1$  with parameters  $\alpha$  and  $\beta$  describing variability in infectivity A simplification of the exact Beta Poisson model when  $\alpha>1$  and  $\beta>10$   $\alpha$ . In that case, the calculation using  ${}_1F_1$  may be very slow, whereas the simplified model is fast and justified (Teunis and Havelaar, 2000).

For the four default index pathogens, the exact Beta Poisson model is applied using pairs of different  $\mathfrak a$  and  $\mathfrak \beta$  values that reflect variability and uncertainty of infectivity of the pathogen. These parameters sets are included in the program code of QMRAspot as Monte Carlo samples of the dose response parameters a and  $\beta$ , and are included in the code of the tool in a packed form to save memory space.

In the run screen, an "other" pathogen can be selected. Norovirus can be selected from the drop-down list, whereby the built-in parameter values for a and  $\beta$  of 0.044 and 0.50 (Teunis et al., 2008) are used. For any other pathogen (which may also be one of the four default index pathogens), a dose response model can be selected and corresponding (fixed) parameter values can be set (see also section 10).

## 6 Help




Figure 6.1 Help-Messages, top of screen

Under the Help-tab, information is provided about the version of QMRAspot and its version history. Also, brief guidance is provided on how to use QMRAspot. Finally, under Help-Messages, (error) messages generated during a run of QMRAspot are listed. The messages also include the Gamma distribution parameter values of the source water concentrations and of all influent and effluent concentrations.

The (error) messages or warnings that may occur are explained here below.

#### {FindMinimum::Istol}

The command FindMinimum may produce the message Istol:

The objective function does not have a smooth minimum.

The line search decreased the step size to within tolerance specified by AccuracyGoal and PrecisionGoal but was unable to find a sufficient decrease in the function. More than MachinePrecision digits of working precision may be needed to meet these tolerances.

The algorithm for finding the maximum likelihood (minimum deviance) parameters found an optimum, but could not establish the precision of the estimated parameter values (because the target function had a too small gradient). Care should be taken to check the appropriateness of the resulting distributions.

#### {FindMinimum::sdprec}

The command FindMinimum may produce the message sdprec: Line search unable to find a sufficient decrease in the function value with MachinePrecision digit precision. This is a similar warning to the previous one. {FindMinimum::nrnum}

The command FindMinimum may produce the message nrnum: This usually implies that the starting values did not lead to finding the appropriate minimum. In this the reported distribution parameter values may not be correct.

In all these cases it is strongly recommended to scrutinize the underlying data. For example, one may test whether parts of the data provide acceptable fitting or insights into the cause of the warning messages.

Support in handling error and warning messages can be provided by sending an email to <a href="mailto:QMRAspot@rivm.nl">QMRAspot@rivm.nl</a>.




Figure 6.2 Help-Messages, scrolled down

Of all raw data, the number of data with counts between 0 and 100, 100 and 200, 200 and 300, 300 and 1000 and more than 1000 are listed. This is done to highlight unrealistically high counts in the data.

## 7 Fitting distributions to the data

#### 7.1 General

All raw data sets should include three or more samples: smaller data sets are ignored and parameters are not estimated. Counts in QMRAdata.xls may only be integers, any non-integer counts are rounded. QMRAspot does not generate a message that it has rounded non-integers.

Distributions are fitted to the data minimizing deviance functions to obtain a maximum likelihood estimate. This procedure finds optimum parameter values by maximizing the likelihood (or posterior probability) of the selected model with the observed data. Maximum Likelihood Estimation (MLE) is a standard approach to parameter estimation and inference in statistics. MLE has many optimal properties in estimation: sufficiency (complete information about the parameter of interest contained in its MLE estimator); consistency (true parameter value that generated the data recovered asymptotically, i.e. for data of sufficiently large samples); efficiency (lowest-possible variance of parameter estimates achieved asymptotically); and parameterization invariance (same MLE solution obtained independent of the parameterization used) (Myung, 2002).

#### 7.2 Source water concentration

Monitoring of the surface water should be aimed at achieving a representative quantification of the numbers of pathogenic microorganisms in the source water, considering seasonal variability, as well as short term fluctuations of pathogen concentrations (Westrell et al., 2006).

If microbial particles in water are homogeneously distributed, then the counts n within each sample of size V are Poisson distributed. Because the concentration is assumed to be gamma distributed among samples, the counts  $n_1...n_N$  of N samples with samples sizes  $V_1...V_N$  have a Negative Binomial distribution with parameters r and  $1/(1+\lambda V_i)$  (Teunis et al., 2009).

A Gamma distribution is used to describe the variability between concentrations in samples taken at different times. A Gamma distribution is a distribution that arises naturally in processes for which the waiting times between events are relevant. It can be thought of as a waiting time between Poisson distributed events. A Gamma distribution has a mean value equal to  $r\lambda$  and a variance equal to  $r\lambda^2$ . Parameter r is a scale parameter and  $\lambda$  a shape parameter.

A Negative Binomial distribution is a discrete probability distribution of the number of successes in a sequence of binomial trials before a specified (non-random) number of failures (r) occurs. For example, if we define a "1" as failure, and all non "1"s as successes, and we throw a dice repeatedly until the third time "1" appears (r = three failures), then the probability distribution of the number of non-"1"s that had appeared will be negative binomial (http://en.wikipedia.org/wiki/Negative binomial distribution).

To summarize, the choice for a Gamma distribution follows from the concept that concentration varies, that it is related to the negative binomial distribution, and that the maximum likelihood estimation is simple.

Parameters r and  $\lambda$  are estimated by minimizing the following deviance function:

$$L[r,1/(1+\lambda V_i)] = -2\ln \prod_{i=1}^{n} [f(r) \operatorname{NegBin}(n_i, V_i | r, s/(s+V_i \lambda))]$$
(1)

where  $f(r)=1-\phi(x-20)/20$  is a prior function for the shape factor r, with  $\phi(x-20)/20$  the cumulative normal distribution with a mean and variance of 20. This function can be interpreted as very flat, prior to that preventing extremely small values of r from occurring, thereby facilitating robust parameter estimation without strongly affecting the estimates. s is a scaling factor. If the mean sample volume is less than one litrr then s=0.001. This scaling avoids computational underflows.

In QMRAdata.xls presence/absence data may be given for any microorganism, although this is usually only the case for *Campylobacter*. These observations are used to calculate a concentration for each actual taken sample by minimizing the following deviance function:

$$L[c,V] = -2\ln \prod_{i=1}^{n} \binom{n_i = 0 \Rightarrow Pois[0|V_i]}{n_i > 0 \Rightarrow 1-Pois[0|V_i]}$$
(2)

where *Pois* denotes Poisson distribution, c is concentration and  $V_i$  is the sample size of the i-th sample.

Subsequently, a Gamma distribution with parameters r and  $\lambda$  is fitted to the concentration data by minimizing the following deviance function:

$$L[r,\lambda] = -2\ln \prod_{i=1}^{n} \begin{pmatrix} c_i = 0 \Rightarrow & \text{GammaCDF}[L[c,V] = \chi^2_{95\%}(\text{df} = 1) | \ln r, \ln \lambda] \\ c_i = \infty \Rightarrow & 1 \text{-GammaCDF}[L[c,V] = \chi^2_{95\%}(\text{df} = 1) | \ln r, \ln \lambda] \\ 0 < c_i < \infty \Rightarrow & \text{Gamma[ln r, ln } \lambda] \end{pmatrix}$$
 (3)

Where  $c_i$  is the concentration of the *i*-th sample. Gammacdf denotes a cumulative function of the Gamma distribution.  $L(c,V) = X^2_{95\%}(df=1)$  is the root of the likelihood function equal to the 95-percentile of a Chi-squared distribution with one degree of freedom (df). Note that in equation (3) the zero and infinity sign were exchanged compared to Schijven et al. (2011). This has been corrected in QMRAspot as of version 1.2.

If raw data of the index pathogen in the source water consist of non-detects only, a Gamma distribution with parameters r=0.01 and  $\lambda=1/(\Sigma V_i+0.01)$  is assumed.

## 7.3 Recovery, R

In order to estimate the recovery efficiency of the detection method of the index pathogen, samples are spiked with a known number of the specified indicator organisms. After processing the samples, a fraction of the spiked organisms will be recovered. The data on the initial spike and on the recovery are paired per experiment. The recovered fraction is assumed to be Beta distributed with parameters a and  $\beta$ . Estimation of these parameters by means of the paired Beta model is explained in detail by Teunis et al. (1999, 2005, 2009). If recovery data are lacking, then QMRAspot applies R=1. Potentially this leads to

underestimation of risk and this should be addressed in the discussion of the OMRA results.

A beta distribution is a family of continuous probability distributions defined on the interval [0, 1] and parameterized by two positive shape parameters, denoted by a and  $\beta$ , that appear as exponents of the random variable and control the shape of the distribution

(http://en.wikipedia.org/wiki/Beta\_distribution). The interval [0,1] applies because recovery is a fraction, i.e. its value lies between 0 and 1.

#### 7.4 Treatment, z

Here, we assume that treatment z is in effect, implying that microorganisms are removed and thus  $0 \le z \le 1$ . It is assumed that microorganisms passing treatment do so independently with a probability or fraction z. This may be modelled as a binomial process, either with paired or unpaired samples (Teunis et al., 1999, 2009), or as the ratio of the Gamma distributed effluent concentrations / the Gamma distributed influent concentrations. Collection of paired data from a treatment step requires exact timing of the sampling. The pairing may be lost if mixing occurs during treatment. Residence times in treatment may vary from a few hours to several days. In many cases, even with short residence times and samples of influent and effluent collected on the same day, pairing is not evident

In QMRAspot under Parameter settings-Treatment, a treatment step can be set to use the raw data from a spreadsheet or to enter characteristic parameter values, or to switch off the treatment step. Depending on these settings and on the raw data values, treatment fraction z may be equal to one, be described by a Beta distribution or, as the ratio of the Gamma distributed concentrations of influent and effluent. This is a so-called type II Beta distribution or F-distribution (Teunis et al., 2009).

Table 1 lists all the possible ways of how z is calculated (model selection) and what output is produced.

Commonly, drinking water companies prefer to characterize each treatment step separately.

However, treatment steps may be combined. In fact, one could combine all treatment steps into one if desired. For example, treatment step 2 and 3 are combined to step 2+3 by analysing the influent data of step 2 and the effluent data of step 3. This needs to be set as such in the data spreadsheet (section 2.5). Of course, influent and effluent data need be associated with each other. A reason for combing treatment may be that the effluent data of step 2 consist only of non-detects, whereas, due to taking larger samples, in the effluent sample of step 3, microorganisms were detected.

## 1 Table 7.1 Model selection for treatment

Setting	Data or condition	Calculation (for explanation of Monte Carlo simulation, see section 8)	Output
Log <sub>10</sub> z,α or α,β	No data needed	$\beta = \alpha/10 ^{\rm log_{10}z}$ -a. a and $\beta$ are the parameter of the Beta distributed treatment fraction z. With RandomReal[BetaDistribution[a, $\beta$ ], 10 000] ten thousand Beta distributed fractions are generated.	Histogram of Beta distributed fractions. The associated histogram heading includes "Parameter values set" and "Unpaired beta model"
Off	No data needed	Treatment fraction z is set equal to one, implying no pathogen removal by this treatment step.	No histogram. Message: "Switched off"
Raw data	No data in the spreadsheet	Treatment fraction z is set equal to one, implying no pathogen removal by this treatment step.	No histogram. Message: "No data"
Raw data	Less than three influent or effluent raw data	Treatment fraction z is set equal to one, implying no pathogen removal by this treatment step.	No histogram. Message: "Too few data"
Raw data	Only non- detects in the influent and effluent data, so zero counts in all samples	Treatment fraction z is set equal to one, implying no pathogen removal by this treatment step.	No histogram. Message: "Too uncertain for use in QMRA "
Raw data	Only non- detects in the influent, but detection in effluent data	Treatment fraction z is set equal to one, implying no pathogen removal by this treatment step. Note that if the sample size of influent and effluent samples are similar, this could imply that z>1.	No histogram. Message: "Too uncertain for use in QMRA"
Raw data	Detection in influent samples and parameter r of the Gamma distributed influent concentration ≥5	If parameter r of the Gamma distributed influent concentration $\geq 5$ , implying near constant influent concentration, then the unpaired beta model has difficulty in finding a solution. To circumvent this technical inconvenience, the gamma ratio model is applied to generate ten thousand fractions. Gamma distribution parameters r (shape) and $\lambda$ (scale) of the influent (i) and effluent (e) concentrations are estimated in the same as for the pathogen concentration in the source. Then z= RandomReal[GammaDistribution[re, $\lambda_e$ ], 10 000] / RandomReal[GammaDistribution[ri, $\lambda_i$ ], 10 000]	Histogram of Ratio of Gamma distributed influent and effluent concentrations. The associated histogram heading includes "High uncertainty".
Raw data	Detection in influent samples, not necessarily in effluent samples	Estimation of the parameters α and β by means of the unpaired Beta model is explained in detail by Teunis et al. (1999, 2009).  With RandomReal[BetaDistribution[α, β], 10 000] ten thousand Beta distributed fractions are generated.	Histogram of Beta distributed fractions. The associated histogram heading includes "Unpaired beta model".

#### 8 Monte Carlo simulation

In QMRAspot, Monte Carlo simulation is used to generate data from probability distributions that describe pathogen concentrations in the source water, the fractions of microorganisms that pass drinking water treatments, and consumption of unboiled drinking water per person per day. The generated random data are used to calculate exposure and infection risk.

The following simple example with Mathematica code illustrates how such Monte Carlo simulation works.

Consider ten throws with a dice, the values are stored in vector a:

a=RandomInteger[{1,6},10]

{1,2,3,1,4,5,1,1,2,5}

In the same way vector b contains the values of ten other throws:

b=RandomInteger[{1,6},10]

{2,5,5,6,3,1,4,1,6,1}

Vector ab is the product of a and b, in which the *i*-th element of a is multiplied with the corresponding *i*-th element of b:

a b

{2,10,15,6,12,5,4,1,12,5}

In QMRAspot, 10,000 Monte Carlo (MC) samples are generated from all distributions. This number of Monte Carlo samples is sufficient (Teunis and Havelaar, 1999). The source water concentration of the index pathogens,  $C_{source}$ , is Gamma distributed with parameters r and  $\lambda/s$ . For recovery R, and treatment steps  $z_1...z_6$ , Beta-distributed MC samples are generated.

For the unboiled drinking water consumption, W litre, MC samples of lognormal distributions are generated, but in case of the WHO-data set, a fixed value of 2 litres is used. Monte Carlo samples of the dose response parameters a and  $\beta$  are provided as pregenerated data and included in the tool in a packed form to save memory space.

## 9 Exposure and infection risk

Exposure to the index pathogens is a given as the dose  $D=C_{source}W$ , the number of ingested index pathogens per person per day and is calculated by multiplying the MC data of the source, recovery and treatment data (maximum 6 treatment steps):

$$D=C_{\text{source}} \frac{1}{R} \prod_{i-1}^{6} z_i W$$
(4)

Infection risk per person per day is calculated by applying the hypergeometric dose-response relation with Beta-distributed dose response parameters a and  $\beta$  (Teunis and Havelaar 2000). The formula for calculating the risk of infection for a specific dose is calculated as follows (Teunis and Havelaar 2000):

$$P_{\text{inf,person,day}} = 1_{-1} F_1(\alpha, \alpha + \beta; -D)$$
 (5)

where  ${}_{1}F_{1}$  is the confluent hypergeometric function.

The probability that a person acquires no infection on the i-th day equals

The probability that a person acquires no infection on any day in a year equals

$$(1-P_{\text{inf,person,day,1}})(1-P_{\text{inf,person,day,2}})...(1-P_{\text{inf,person,day,365}}) = \prod_{i=1}^{365} (1-P_{\text{inf,person,day,i}})$$
(8)

Note that the daily risks are assumed to be independent of each other. Thus, the distribution of the probability or risk that a person per year is infected at least once, denoted simply as the infection risk per person per year, is calculated from MC sampling from the daily infection risk (Teunis et al., 1997):

$$P_{\text{inf,person,year}} = 1 - \prod_{i=1}^{365} \left( 1 - P_{\text{inf,person,day,i}} \right)$$

$$\tag{9}$$

## 10 Dose response data

Applied dose response relations were generally derived from studies in which a specific strain of the index pathogen was given to human volunteers (Teunis et al., 1996, 2002a, 2002b). However, one pathogen strain does not represent the suite of strains that may occur in source waters for drinking water production. A hierarchical dose response relation, as was performed for multiple isolates of *Cryptosporidium parvum*, produced estimates that differed greatly between isolates (Teunis et al. 2002a). Predictions based on multilevel dose response relations may aid probabilistic risk assessments such as those presented here to properly reflect the variation in pathogen strains. Moreover, dose response data from outbreaks may inform the dose response relation as was shown for *Campylobacter* (Teunis et al. 2005b, 2008a, 2010, Thebault et al. 2013). Such additions, both hierarchical analysis and the use of outbreak data, could aid the estimation of the enterovirus dose response relation for which now the rotavirus dose response relation is used. Data are currently not available for this type of analysis, and additional research is required.

A brief overview of published dose response information, including some statistics) is given in Table 10.1. Dose response assessments based on multilevel models produce estimates for the dose response parameters ( $\alpha$ ,  $\beta$ ) as (joint) distributions. For these analyses, no alpha or beta estimates are given in Table 10.1. The latter applies to the four index pathogens: enterovirus, *Campylobacter, Cryptosporidium* and *Giardia*.

2 Table 10.1 Dose response data

Name	alpha	beta	Low dose	ID50	Ref	
			inf			
			= a/(a+β)			
Vibrio cholerae	0.508	7.52x10 <sup>7</sup>	7.10x10 <sup>-9</sup>	2.13x10 <sup>8</sup>	Teunis et al. (1996)	
	0.164	0.136	3.79x10 <sup>-2</sup>	$1.16 \times 10^{2}$	Teunis et al. (1996)	
Salmonella Typh/Ent	8.53x10 <sup>-3</sup>	3.14	2.71x10 <sup>-3</sup>	6.65	Teunis et al. (2010)	
Campylobacter jejuni	0.024	0.011	0.685	1.29	Teunis et al. (2005)	
E.coli O157:H7	-	-	6.80x10 <sup>-3</sup>	6.18x10 <sup>2</sup>	Teunis et al. (2008a)	
Cryptosporidum parvum	-	-	0.213	17.6	Teunis et al. (2002)	
Cryptosporidium hominis	8.37x10 <sup>-11</sup>	2.62x10 <sup>-11</sup>	0.762	1.07	Chappel et al.( 2006)	
Giardia lamblia	r=1.9	9x10 <sup>-2</sup>	1.99x10 <sup>-2</sup>	34.83	Teunis et al. (1996)	
Entamoeba coli	0.106	0.295	0.193	2.36x10 <sup>2</sup>	Teunis et al. (1996)	
Norwalk virus (Se+)	0.040	0.055	0.421	18.0	Teunis et al. (2008b)	
Norovirus (Se+)	-	-	0.515	16.2	Thebault et al. (2013)	
Rotavirus	0.167	0.191	0.466	3.46	Teunis et al. (2000)	
Echovirus	0.401	227.2	1.77x10 <sup>-3</sup>	$1.05 \times 10^{3}$	Teunis et al. (1996)	

The parameters for Norwalk virus (Se+) are built in the code of QMRAspot.

Low dose: If the exposure or dose is far less than one,  $1_{-1}F_1(a,a+\beta;-D)$  can be simplified to  $(a/(a+\beta))D$ .

## 11 Example peak concentration in source water

The following example shows the effect of a peak value of the pathogen concentration in the source water. Fig 11.1 shows the values (top left) and distribution (bottom left) of the concentration of Enterovirus in the source water taken from the reference QMRA-spreadsheet. Of the 35 samples of approximately 500 litres, five were found positive. In each of these five samples, only one virus particle was detected. All five concentrations are well above the 95-percentile, and should, therefore, be scrutinized whether they are peak concentrations (outliers). The right panel shows the same data again, but sample 35 was changed from a non-detect to a sample in which ten virus particles were counted. In this case, all the samples with only one virus particle now fall within the 95%-area and the only outlier is the simulated peak concentration.

Note how the Gamma distribution parameter values changed between the two cases. The distribution for the case with the one high peak concentration value is wider (Gamma parameter r is considerably smaller) and the mean concentration value has increased by  $0.5 \log_{10}$  (about three times higher;  $\lambda$  is almost five times higher). Obviously, this implies that the infection risk will also be about three times higher. This example shows that missing a peak concentration may lead to a significant underestimation of the risk, of which one is unaware.




Figure 11.1 Effect of a peak concentration in the source water. Top left: Enterovirus concentrations in 35 samples. Top right: One non-detect replaced by a ten times higher peak concentration. Bottom left and right: corresponding Gamma distributions.

## 12 Example treatment

#### 12.1 General

The following example is a QMRA with only one positive sample in the effluent of the last treatment step. These are the data from a drinking water company. Only the data and analysis of the last treatment step z2 and the risk outcome are shown. This QMRA is compared with the case in which the effluent of the last treatment step only contained non-detects.

This situation is encountered regularly. The percentage of non-detects obviously is the highest in the last treatment step of a drinking water production.

#### 12.2 Treatment z2, one positive sample in effluent




Figure 12.1 Data and analysis of treatment z2 with one positive effluent sample

Over a period of six years, more than 1200 samples were analysed in the influent and effluent. Only one sample in the effluent was found positive. Sample size of the effluent samples was ten times that of the influent samples. The parameter values estimated in QMRAspot for the Beta distribution that describes the variability of z2 are  $\alpha = 0.011$  and  $\beta = 68$ . The low value of  $\alpha$  reflects the very wide distribution of z2. Although more than 1200 samples were analysed, only one bacteria colony was counted, so there is very limited information to estimate the removal efficiency of the treatment step. In such a wide distribution, the mean value of z carries little information on the reduction of the risk of exposure and infection.

# 12.3 Treatment z2, only non-detects in effluent




Figure 12.2 Data and analysis of treatment z2 with one non-detects in effluent

Having only non-detects in the effluent implies that there is even less information available to estimate removal efficiency of the treatment step, which is reflected by an extreme low value  $\alpha = 0.000011$  and a histogram showing a uniform distribution that spans a huge range of  $\log_{10}z$ . Obviously, the mean value of z is here meaningless.

## 12.4 Infection risk, one positive sample in effluent




Figure 12.3 Infection risk (treatment z2 with one positive effluent sample)

Figure 12.3 shows that in the scenario with only one positive sample (count=1), the 95-percentile is about  $10^{-4}$  per person per year.

#### Run Parameter settings | Enterovirus | Campylobacter | Cryptosporidium | Giardia | Report | Help Data scheme | Source water | z1 | z2 | z3 | z4 | z5 | z6 | Exposure | Infection risk | Summary table | Sensitivity analysis Campylobacter-Infection risk Monte Carlo histograms 0.20 logx -7.2 0.5 0.15 P 0.10 0.0 5% Median 1edian -0.5 0.05 -1.0 0.00 Log<sub>10</sub> Infection risk/person/day Log<sub>10</sub>Infection risk/person/year Campylobacter -10-15 -20 Campylobacter

## 12.5 Infection risk, only non-detects in effluent

Figure 12.4 Infection risk (treatment z2 with one non-detects in effluent)

Figure 12.4 shows that in this scenario with only no-detects in the effluent of treatment z2, the 95-percentile is slightly higher compared to the other scenario, namely  $10^{-3.6}$  per person per year. Most striking is the extremely wide range of the predicted infection risk. Compared to the alternative scenario, the minimum risk is much lower. The maximum risk is slightly higher than in the other scenario. These differences are all due to the fact that in the scenario with only non-detects, estimates are much more uncertain.

Using the data of z1 in the example, a number of scenarios were followed in which the number of detects in the effluent of z1 was varied to investigate the changes in the parameter values of the Beta distribution that describes z1, and to make a comparison with the point estimate of z1.

#### 12.6 Monte Carlo versus Point estimates

In understanding the outcomes of the QMRA for drinking water, one is inclined to make a comparison with the calculations based on the mean values of each individual step in the QMRA on the  $\log_{10}$  scale, for example, as in calculating total treatment. The above example is used as an example to compare the point estimates with the mean value of the Monte Carlo data for treatment step z2 (Table 12.1).

The point estimate of fraction z2 is about  $2\log_{10}$  lower (more removal) than the mean of the Monte Carlo data that are sampled from the corresponding Beta distribution. Recall that for z2, 1248 influent and 1248 effluent samples were analysed. In only one effluent sample, one bacterium was counted. It shows that for three cases with one positive sample in which 1, 2 or 3 bacteria were counted respectively, that parameter  $\alpha$  remains the same, but parameter  $\beta$  gradually decreases (Table 12.1). The mean of z (and of  $\log_{10}z$ ) remains approximately constant, as well as the shape of the Beta distribution, whereas the point estimates linearly increases with the counted bacteria. The last two scenarios demonstrate the importance of the distribution. If two or three counts are distributed over two and three samples respectively, then the mean of z decreases, whereas the point estimate is unaffected. This demonstrates the superiority of using a distribution to describe the data over a point estimate.

## 3 Table 12.1 Comparison of Monte Carlo data with the point estimate of fraction z for treatment step z2.

Scenario for effluent of z2	α	ß	Mean of z	Mean of	Point estimate	log <sub>10</sub> of point estimate
	0.00007	0.10	0.00070	log <sub>10</sub> z		
Only non-detects	0.000097	0.18	0.00079	-3.1	0	-∞
1 count	0.011	68	0.00016	-3.8	1.7x10 <sup>-6</sup>	-5.8
2 counts/one sample	0.011	65	0.00017	-3.8	3.4x10 <sup>-6</sup>	-5.5
3 counts/one sample	0.011	63	0.00017	-3.8	5.1x10 <sup>-6</sup>	-5.3
2 samples with one count	0.0046	11	0.00046	-3.3	3.4x10 <sup>-6</sup>	-5.5
3 samples with one count	0.0062	8.1	0.00071	-3.1	5.1x10 <sup>-6</sup>	-5.3

## 13 Contact

For user questions or feedback, please contact <a href="mailto:QMRAspot@rivm.nl">QMRAspot@rivm.nl</a> or <a href="mailto:Jack.Schijven@rivm.nl">Jack.Schijven@rivm.nl</a>.

To obtain a dataset with the pairs of a and ß parameter values for the dose response of the index pathogens as built in the code of QMRAspot, a request should be sent to <a href="mailto:Peter.Teunis@rivm.nl">Peter.Teunis@rivm.nl</a>.

Discussions, feedback, questions and answers can be also found in the LinkedIn group op QMRAspot.

#### 14 References

- Anonymous. (2001) Besluit van 9 januari 2001 tot wijziging van het waterleidingbesluit in verband met de richtlijn betreffende de kwaliteit van voor menselijke consumptie bestemd water. Staatsblad van het Koninkrijk der Nederlanden, 31:1-53.
- Anonymous. (2005) VROM-Inspectierichtlijn Analyse microbiologische veiligheid drinkwater, VROM-Inspectie, Artikelcode 5318.
- Benford, D. (2001). Principals of Risk Assessment of Food and Drinking Water Related to Human Health. ILSI Europe concise monograph SERIES International Life Science Institute, Belgium pp.1-43.
- Chappell CL, Okhuysen PC, Langer-Curry R, Widmer G, Akiyoshi DE, Tanriverdi S, Tzipori S (2006). *Cryptosporidium hominis*: experimental challenge of healthy adults. American Journal of Tropical Medicine and Hygiene 75(5):851-857.
- De Roda Husman, A. M., Lodder, W.J., Rutjes, S.A., Schijven, J.F. and Teunis, P.F.M. (2009) Long-term inactivation study of three enteroviruses in artificial surface and groundwaters, using PCR and cell culture. Applied Environmental Microbiology 75(4):1050–1057.
- Haas, C.N. "How to average microbial densities to characterize risk." Water Research 30.4 (1996): 1036-1038.
- Haas, C.N. and Eisenberg, J.N.S. (2001). Risk Assessment. In Fewtrell, L. and Bartram, J. (ed.), Water Quality: Guidelines, Standards and Health. IWA Publishing, London, p. 161-183.
- Haas, C.N., Rose, J.B. and Gerba C.P. (1999). Quantitative microbial risk assessment. Wiley, New York, USA.
- Lodder, W.J., van den Berg, H.H., Rutjes, S.A. and de Roda Husman, A.M. (2010) Presence of enteric viruses in source waters for drinking water production in The Netherlands. Applied Environmental Microbiology 76(17):5965-71.
- Medema, G.J., Hoogenboezem, W., van der Veer A.J., Ketelaars, H.A.M., Hijnen, W.A.M. and Nobel, P.J. (2003) Quantitative Risk Assessment of *Cryptosporidium* in surface water treatment. Water Science & Technology 47:241-247.
- Myung IJ. Tutorial on maximum likelihood estimation. J. Math Psych, 2003:47,90–100.
- Rutjes, S.A. and de Roda Husman, A.M. (2004) Procedure for virus detection in water for use in the Dutch drinking-water guidelines. RIVM Report 330000007, 38 p in Dutch.
- Rutjes, S.A., Lodder, W.J., van Leeuwen, A.D. and de Roda Husman, A.M. (2009) Detection of infectious rotavirus in naturally contaminated source waters for drinking water production. Journal of Applied Microbiology 107(1):97-105.
- Schets, F.M., Medema, G.J. and Schijven, J.F. (2004) The recovery of the detection method for *Cryptosporidium* and *Giardia* in water. RIVM Report 330000008, 71 p in Dutch.
- Schets, F.M., van Wijnen, J.H., Schijven, J.F., Schoon, H. and de Roda Husman, A.M. (2008) Monitoring of waterborne pathogens in surface waters in amsterdam, the Netherlands, and the potential health risk associated with exposure to *Cryptosporidium* and *Giardia* in these waters. Applied Environmental Microbiology 74(7):2069-78.

- Schijven, J. F., Teunis, P. F., Rutjes, S. A., Bouwknegt, M., & de Roda Husman, A. M. (2011). QMRAspot: A tool for Quantitative Microbial Risk Assessment from surface water to potable water. *Water research*, 45(17), 5564-5576.
- Teunis PFM, Kasuga F, Fazil A, Ogden ID, Rotariu O, Strachan NJC (2010). Doseresponse modeling of *Salmonella* using outbreak data. International Journal of Food Microbiology 144(2):243-249.
- Teunis PFM, Moe CL, Liu P, Miller SE, Lindesmith L, Baric RS, Le Pendu J, Calderon RL (2008b). Norwalk virus: how infectious is it? Journal of Medical Virology 80(8):1468-1476.
- Teunis PFM, Ogden ID, Strachan NJC (2008a). Hierarchical dose response of *E. coli* O157:H7 from human outbreaks incorporating heterogeneity in exposure. Epidemiology and Infection 136(6):761-770.
- Teunis PFM, van der Heijden OG, van der Giessen JWB, Havelaar AH (1996). The dose response relation in human volunteers for gastro-intestinal pathogens. RIVM report 284550002.
- Teunis, P.F.M, Chappell, C.L. and Okhuysen, P.C. (2002a) *Cryptosporidium* dose response studies: variation between isolates. Risk Analysis 22(1):175-183.
- Teunis, P.F.M, Rutjes, S.A., Westrell, T. and de Roda Husman, A.M. (2009) Characterization of drinking water treatment for virus risk assessment. Water Research 43(2):395-404.
- Teunis, P.F.M. (1997) Infectious gastroenteritis opportunities for dose response modelling. RIVM report 284550003.
- Teunis, P.F.M. and Havelaar, A.H. (2000) The Beta Poisson model is not a single hit model. Risk Analysis 20:511-518.
- Teunis, P.F.M. and Havelaar, A.H. (1997) Cryptosporidium in drinking water: evaluation of the ILSI quantitative risk assessment framework.RIVM report 28455006.Teunis, P.F.M., Chappell, C.L. and Okhuysen, P.C. (2002b) Cryptosporidium dose-response studies: variation between hosts. Risk Analysis 22(3):475-485.
- Teunis, P.F.M., Evers, E.G. and Slob, W. (1999) Analysis of variable fractions resulting from microbial counts. Quantitative Microbiology 1:63-88.
- Teunis, P.F.M., Lodder, W.J., Heisterkamp, S.H. and de Roda Husman, A.M. (2005a) Mixed plaques: statistical evidence how plaque assays may underestimate virus concentrations. Water Research 39(17):4240-50.
- Teunis, P.F.M., Medema, G.J., Kruidenier, L. and Havelaar, A.H. (1997). Assessment of the risk of infection by *Cryptosporidium* or *Giardia* in drinking water from a surface water source. Water Research 31(6):1333-1346.
- Teunis, P.F.M., Van den Brandhof, W., Nauta, M., Wagenaar, J., Van den Kerkhof, H. and Van Pelt, W. (2005b). A reconsideration of the *Campylobacter* dose-response relation. Epidemiology and Infections 133(4):583-92.
- Teunis, P.F.M., Van der Heijden, O.G., Van de Giessen, J.W.B and Havelaar, A.H. (1996) The dose-response relation in human volunteers for gastro-intestinal pathogens. RIVM Report 284550002.
- Teunis, P. F. M., and Havelaar, A. H. (1999). Cryptosporidium in drinking water: evaluation of the ILSI quantitative risk assessment framework. RIVM Report 284550006.
- Thebault A, Teunis PFM, Le Pendu J, Le Guyader F, Denis J-B (2013). Infectivity of GI and GII noroviruses from oyster related outbreaks. Epidemics 5(2): 98-110.
- USEPA (2006) Economic Analysis for the Final Ground Water Rule. United States Environmental Protection Agency. 2006; EPA 815-R-06-014. <a href="https://www.epa.gov/safewater">www.epa.gov/safewater</a>
- Westrell, T., Teunis, P., Van den Berg, H., Lodder, W., Ketelaars, H., Stenström, T.A. and de Roda Husman A.M. (2006) Short- and long-term variations of

norovirus concentrations in the Meuse river during a 2-year study period. Water Research 40: 2613-2620.

WHO 2004. Guidelines for Drinking-water Quality. Third edition, World Health Organisation, Geneve, Zwitserland. ISBN 9789241546386

WHO 2011. Guidelines for Drinking-water Quality. Fourth edition, World Health Organisation, Geneve, Zwitserland. ISBN 978 92 4 154815 1.